

ELECTRIC TOOTH BRUSH (on batteries)

INSTRUCTION MANUAL



Dear consumer, Please read these instructions carefully before use in order to avoid damage due to incorrect use and to better understand the functions of this product. Please pay particular attention to the safety guide. If you pass the device on to a third party, these Operating Instructions must also be handed over. Please keep them for future references.

IMPORTANT SAFEGUARDS

- ❑ Please read these instructions carefully before use.
- ❑ Keep present instructions for future references.
- ❑ Only use the appliance for domestic purposes and household or similar applications such as:
 - staff kitchen areas, in shops, offices and other working environments;
 - By clients in hotels, motels and other residential type environments;
 - Environment type like :
 - Farm Houses • Bed and Breakfast.
- ❑ Only use the appliance in the way indicated in these instructions.
- ❑ Never use this unit near bathtub, shower, wash-hand basins or other containers with water.
- ❑ Never use this apparatus near water projections.
- ❑ Never use this apparatus with the wet hands.
- ❑ Inform potential users of these instructions.
- ❑ Never leave the appliance unsupervised when in use.
- ❑ The appliance must be used only for its intended purpose. No responsibility can be taken for any possible damage caused by incorrect use or improper handling.
- ❑ Children and persons with disabilities are not always aware of the danger posed by the use of appliances.
- ❑ This device can be used by children aged 8 years and older and adults, with reduced physical, sensory or mental capabilities or lack of experience and knowledge, if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- ❑ This device is not a toy, do not leave small children or people with disabilities to play with.
- ❑ Cleaning and maintenance must not be made by children without supervision.

- In order to ensure your children's safety, please keep all packaging (plastic bags, boxes, polystyrene etc.) out of their reach.
- Caution ! Do not allow small children to play with the foil :
THERE IS A DANGER OF SUFFOCATION !
- Never immerse the appliance in water or any other liquid for any reason whatsoever.
- Never place it into the dishwasher.
- Never use the appliance near hot surfaces.
- Never use the appliance outside and always place it in a dry environment.
- Never use accessories, which are not recommended by the producer. They could constitute a danger to the user and risk to damage the appliance.
- Incorrect operation and improper use can damage the appliance and cause injury to the user.
- This device, designed for home use, complies with standards relating to this type of product.
- The item is not to be used, if it has been dropped or if there are visible signs of damage.
- Regarding the instructions for cleaning, thanks to refer to the below paragraph of the manual.
- This appliance is only intended for household use.

USE

Put the toothpaste on the brush and gently insert it into the mouth.

Turn the switch on to start cleaning the teeth and tongue.

NOTE: During the first use, it is possible to observe a slight bleeding of the gums. It should stop after several uses.

If bleeding persists after two weeks of use, contact the dentist.

After use, rinse your mouth and rinse thoroughly the brush head, remove excess water, dry and store the brush.

It is recommended : to replace the head at least every three months.

The monitoring of an adult if the brush is used by a child
(not suitable for children under 3 years).

The use of the electric brush by children must always be accepted beforehand by a dentist or a paediatrician.

CLEANING & MAINTENANCE

After rinsing the toothbrush, remove the excess water and dry the brush with a cloth or towel and store the brush.

Never use abrasive products.

Make sure the brush handle is completely dry before replacing the batteries.

SPECIFICATIONS

Works with 2 x 1.5V AA or LR6 batteries

Made in PRC

Features may change without prior notice

This apparatus satisfied EC directives, it was controlled according to all the European current directives.

This apparatus was designed and manufactured in respect of the last safety regulations and technical specifications.

Imported by Sotech International

" LE PERIPOLE " N°A-107

**33, Avenue du Marechal de Lattre de Tassigny
94120 Fontenay sous bois - France**